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GB Instruction Manual

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BRESSER Science MTL-201
Metallurgical Reflection-type Microscope

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CAUTION!

To work with this microscope, sharp and pointed aids are being used. Please take care that this microscope and its accessories are stored at a place out of reach of children. Let children only work with this microscope under an adult's supervision! Keep packing material (plastic bags etc.) away from children!

Note!

If you have any complaints or queries please first contact your national service centre by telephone. The address is included in these instructions.

Only for EU countries

Do not dispose of electric equipment together with household waste material!

In observance of European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and its implementation in accordance with national law, electric equipment that have reached the end of their life must be collected separately and returned to an environmentally compatible recycling facility.

Discharged batteries and damaged re-chargable batteries must be disposed of at special battery collection points. Information is available from your local disposal agent or local authority regarding the disposal of devices or batteries manufactured after the 01.06.2006

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Kaufdatum: _____

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Notes

Lined writing area consisting of two columns of horizontal lines.

General Information

About this Instruction Manual

Please read the safety instructions in this manual carefully. To avoid damage to the unit and the risk of injury, only use this product as described in the manual.

Keep the instruction manual handy so that you can easily look up information on all the functions.

DANGER!

You will find this symbol in front of every section of text which deals with the risk of severe injury or even death in the event of improper use.

CAREFUL!

You will find this symbol before every section of text that deals with the risk of minor to severe injuries resulting from improper use.

NOTE!

You will find this symbol in front of every section of text which deals with the risk of damage to property or the environment.

Intended Use

This product is intended only for private use. It was developed for the magnified display of things in nature.

General Warning

DANGER!

Tools with sharp edges and points are often used when working with this device. For this reason, store this device and all accessories and tools in a location that is out of the reach of children. There is a RISK OF INJURY!

DANGER!

This device contains electronic components which operate via a power source (power supply and/or batteries). Do not leave children unattended when using the device. Only use the device as described in the manual, otherwise you run the RISK of an ELECTRIC SHOCK.

DANGER!

Do not expose the device to high temperatures. Use only the supplied power adapter. Do not short circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short circuit, fires or even explosions.

DANGER!

Never bend, pinch, pull or run over the power and connecting cables or extensions/adapters. Protect the cables from sharp edges and heat. Before operation, check the device, cables and connections for damage. Never use a damaged unit or a unit with damaged power cables, etc. Damaged parts must be exchanged by an authorized service centre immediately.

Position your device so that it can be disconnected from the power supply at any time. The wall socket should always be located near the device and be easily accessible, since the plug on the power cord serves as a disconnecting device for the power supply.

DANGER!

Children should only use the device under adult supervision. Keep packaging material (plastic bags, rubber bands, etc.) out of the reach of children! There is a RISK OF CHOKING!

CAREFUL!

Children must not have access to the included chemicals and liquids. Do not drink the chemicals. Wash hands thoroughly with running water after use. In the event of contact with the eyes or mouth, rinse thoroughly with water. In the event of pain, contact a doctor immediately and take the substances with you.

NOTE!

Do not disassemble the device. In the event of a defect, please contact the Service Centre in your country (see the warranty card).

Do not expose the device to temperatures above 40°C.

NOTES on cleaning

Separate the device from the power supply before cleaning (remove plug).

Only use a dry cloth to clean the exterior of the device.

NOTE!

Do not use any cleaning fluid to avoid damaging the electronics.

Clean the lenses (eyepiece and/or lens) with a soft and lint-free cloth only (e.g. microfibre).

NOTE!

Do not apply excess pressure to the cloth so as to avoid scratching the lenses.

To remove more stubborn dirt, moisten the cleaning cloth with an eyeglass-cleaning solution and wipe the lenses gently.

Protect the device from dust and moisture! After use - in particular in situations of high humidity - let the device acclimatize for a short period of time, so that the residual moisture can dissipate. Remove the dust cover and store it in the included bag.

Disposal

Dispose of the packaging materials properly, according to their type (paper, cardboard, etc). Contact your local waste disposal service or environmental authority for information on the proper disposal.

Do not dispose of electronic devices in the household garbage! As per the Directive 2002/96/EC of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner. Empty old batteries must be disposed of at battery collection points by the consumer. You can find out more information about the disposal of devices or batteries produced after 01.06.2006 from your local waste disposal service or environmental authority.

Please take the current legal regulations into account when disposing of your device. You can get more information on the proper disposal from your local waste disposal service or environmental authority.

WARNING

Please disconnect the plug from mains socket before changing the fuse or opening the lamp house to change the lamp. The power cord provided with the equipment has a grounded plug. Always use the power cord with a properly grounded mains socket.

Do not expose the instrument to high temperatures or humidity. Avoid using the instrument in extremely dusty locations.

Operating temperature: 5° C to 35° C.

Operating humidity: 20% to 80% (at 25° C).

Table of Contents

- I. Components
- II. Technical Data
- III. Installation
- IV. Basic Operation
- V. Maintenance, EEC Conformity Explanation, Warranty

I. COMPONENTS

Fig. 1

- 1 Filter seat
- 2 Field diaphragm centering screw
- 3 Polarizer adjustment pole
- 4 Phototube
- 5 Objective
- 6 Stage
- 7 Limit knob (up stop, pre-focusing lever)
- 8 Coarse focusing knob
- 9 Fine focusing knob
- 10 Quick adjustment rod

Fig. 2

- 11 Eyepiece
- 12 Binocular tube set screw
- 13 Field diaphragm adjustment rod
- 14 Aperture diaphragm adjustment rod
- 15 Power supply socket (mains in)
- 16 Power cord
- 17 Power switch
- 18 Brightness control knob
- 19 Knob for tensional adjustment of focusing
- 20 Transversal adjustment knob
- 21 Longitudinal adjustment knob

III. INSTALLATION

1. Remove all parts from their packing materials and retain the packaging in the event you need to transport the product.
2. Loosen the binocular tube set screw. Install the binocular tube to the optical body by removing its dust cap and inserting the dovetailed flange on the underside of the binocular tube into the optical body. Secure the binocular tube with the screw.

Do not release the optical part from your grasp until it is firmly secured to the frame.

3. Install the adjustment rods and the filter seat to the frame.
4. Install the objectives into the nosepiece.
5. Remove the dust caps from the eyetubes and insert both eyepieces into the eyetubes.
6. Connect the power cord to the mains in of the microscope and to a suitable mains socket.

II. TECHNICAL DATA

Eyepieces

Eyepiece type	Magnification	Field of view	Brennweite	Remarks
Wide field eyepiece (WF)	10x	\varnothing / mm 18	f / mm 25	
WF eyepiece with crosshair	10x	18		optional

Objectives

Objective type	Magnification	Numerical aperture	Remarks
Planachromatic objective (without covering glass usable)	4x	NA 0,1	optional
	5x	0,12	
	10x	0,25	
	20x	0,40	
	40x	0,60	
80x	0,80		

Total Magnification

Objectives:	4x	5x	10x	20x	40x	80x
Total Magnification						
Eyepiece:						
	10x	40x	50x	100x	200x	800x

- Power supply: AC 230 V / 50 Hz.
- Illumination: Halogen lamp (6 V, 20 W), brightness adjustable.
- Colour filters: matted glass, blue, green and yellow.

- Large stage (mechanical desk and manually movable).
- Minimum division of fine focusing scale: 0.0007 mm.
- Adjustment range of interpupillary distance: 53 mm - 75 mm.
- Fuse: 250 V, 0.5 A, F.

Reservation of technical alterations and errors

IV. . BASIC OPERATION

1. Illumination Controls (Fig. 3)


Fig. 3

- 1 Power switch 2 Brightness control knob

The power switch (1) for the illuminator and the brightness control knob (2) are located on the microscope base. The electrical system is fuse-protected and the fuse holder is located on the power supply socket.

Turn on the light with the power switch (1). If the light does not appear, check the brightness control (2) to see if it's on a sufficiently low setting.

2. Focusing Controls (Fig. 4)


Fig. 4

- 1 Coarse focusing knob
2 Limit knob (up stop, pre-focusing lever)
3 Knob for tensional adjustment of focusing
4 Fine focusing knob

1) Focusing adjustment: Focusing adjustment is accomplished by using the large coarse adjustment knobs (1) located comfortably on each side of the frame. Fine adjustment is accomplished using the smaller knobs (4) located on the same focusing shaft. This coaxial arrangement allows for easy, precise adjustment without drift or discomfort.

2) Focus control: Turning either of the coarse focusing knobs (1) will raise or lower the stage. The smallest graduation on the index scale of the fine focusing knob (4) is approximately 0.7 μm of vertical.

3) Focusing tension adjustment: The tension of the coarse focusing knob is adjustable and preset by the factory for ease of use. If you wish to adjust the coarse focusing tension, first locate the ringed tension adjustment knob (3). It is located between the frame and the right coarse

focusing knob. Turning the knob toward the rear of the microscope increases the tension, and toward the front of the microscope loosens it. The tension is too high if you experience physical discomfort.

4) Pre-focusing or focus stop control: The use of this feature will insure that the shorter working distance objectives don't contact the stage or specimen when using the microscope. Also, its use simplifies focusing. After focusing on the specimen with the coarse focusing knob by a low magnification eyepiece, turning up and fastening of the lever (2) will set an upper limit on the coarse focusing movement. After changing specimens or objectives, focusing is easily accomplished by rotating the coarse focusing knob to reach the pre-focused position. Then make fine adjustments with the fine focusing knob. Focusing movement with the fine focusing knob isn't affected by using the pre-focusing lever.

3. Interpupillary Distance and Diopter Adjustments (Fig. 5)

5)

Fig. 5

- 1 Binocular tube 2 Diopter adjustment ring

1) Interpupillary distance adjustment: The proper interpupillary distance, or the distance between the eyepieces, is crucial to the comfort of the user. Adjusting the interpupillary distance is accomplished via a „folding“ action of the binocular tube (1), allowing for quick and easy adjustment.

2) Diopter adjustment: The proper correction for the individual vision is accomplished by the diopter adjustment located at the left eyepiece. Using the objective with 40x magnification, close your left eye and bring an image into focus with your right eye only. Once the image is well focused, observe with the left eye and make fine adjustments with the diopter adjustment ring (2) to correct for your vision.

4. Mechanical Stage Controls (Fig. 6)

The mechanical stage control consists of three parts:

Fig. 6

- 1 For quick adjustment, push the lever (1) for moving towards a position.
- 2 Cross moving knob for transversal adjustment (2).
- 3 Lengthwise moving knob for longitudinal adjustment (3).

5. Vertical Illumination Alignment (Fig. 7)


Fig. 7

- 1 Phototube
- 2 Field diaphragm adjustment rod
- 3 Field diaphragm centering screw
- 4 Filter seat
- 5 Aperture diaphragm adjustment rod
- 6 Focusing lens adjustment lever
- 7 Halogen lamp
- 8 Cover
- 9 Screws for cover
- 10 Lamp socket
- 11 Lamp centering screws

- 1) The aperture diaphragm near the lamp house may be adjusted slightly (with 5) to change the contrast.
- 2) To center the incident or vertical lamp:
 - a) Put a piece of paper and place it on the stage. Take out one objective from the nosepiece and rotate the nosepiece until the objective opening hole is in position over the paper.
 - b) Turn the power on and adjust the brightness to establish sufficient light.
 - c) This paper will allow for focusing the lamp filament on it. If the filament image (Fig. 8, No. 1) is not centered in the overpass hole (Fig. 8, No. 2), reduce the aperture diaphragm (with 5), loosen the screws (9), remove the cover (8) and

adjust the lamp centering screws (11) to obtain an uniform image. Adjust the focusing lens lever (6) until a sharp image of the lamp filament is attained. Refer to Fig. 8 for this step.


Fig. 8

- 1 Filament image
- 2 Overpass hole

- 3) Reinstall the objective again.
- 4) Adjust the focusing control until the image is clear. If the brightness of the field not equal, you may adjust the lamp centering screws (11) and the focusing lens lever (6) slightly to make the brightness of the field equal. Then reinstall the cover.
- 5) Reduce the field diaphragm (with 2). If it is not centered, move it to the center using the field diaphragm centering screws (3). (These screws are located on both sides behind small protecting caps, which have to be removed beforehand.) Adjustment can be done by opening the field diaphragm.
- 6) The filter seat contains a green, a yellow, a blue filter and a matted glass. You may drive anyone to the light path in order to obtain a image of special quality.
- 7) The polarization analyzer is used in conjunction with the polarization filter cube, selected by pushing in the analyzer filter rod (Fig. 9, No. 4). Push the polarizer (Fig. 1, No. 3) to the working position, rotate it while watching the brightness in the field and let the polarizer and analyzer stand orthogonally (minimum brightness). In order to remove the analyzer from light path, push the rod for normal upper light (Fig. 9, No. 3) in.
- 8) The rod (Fig. 9, No. 2) to the upper left of the trinocular head is to direct 100% of the image light to the binocular tube (with eyepieces) or to the phototube (Fig. 9, No. 1). It satisfies high-grape photography.

Fig. 9


- 1 Phototube with camera port
- 2 Rod to direct the image to binocular tube (with eyepieces) or to camera port
- 3 Rod for normal upper light
- 4 Rod for analyzer filter

6. Changing the Lamp (Fig. 7)

1. Turn off and unplug the microscope.
2. Loose the screws (9) and remove the cover (8).
3. The halogen lamp (7) can be removed by grasping the bulb and pulling it firmly from its socket (10). Take care not to touch the lamp with bare fingers, as the lamp will be hot.
4. Insert a new halogen lamp into the same socket. When you are installing the new lamp, be careful not to touch the glass with your fingers. The new lamp should be supplied in a plastic protective envelope. If not, use a tissue or another medium to grasp the bulb. This will prevent contamination from your hand from reducing the lamp's intensity and lifetime.
5. Reinstall the cover to the frame.
6. If necessary, you may need to readjust the lamp as previously described.

7. Changing the Fuse (Fig. 10)

- ① Fuse holder
- ② Power supply socket (mains in)

Switch off the power switch (Fig. 2, No. 17) and pull out the plugs of the power cord (Fig. 2, No. 16) from mains socket and from mains in (2) at the microscope. Open the fuse holder: Grip with a

small screwdriver under the lug of the fuse holder (1) and pull it carefully forwards (Fig. 10, arrow). Take out the fuse holder (1). Change the fuse. Install the fuse holder and plug again. The specification of fuse: 250 V, 0.5 A, F.

V. MAINTENANCE

1. Cleaning the frame and the stage

Disconnect the plug from mains socket before cleaning. Clean the frame and the stage with a soft cloth moistened with a mild detergent solution. Be sure that the instrument is dry before using.

2. Cleaning optical parts

The eyepieces and the objectives of the microscope are coated. They should not be wiped because dry dirt or dust may scratch the coating. It is best to remove the parts to be cleaned from the frame prior to cleaning. Always blow loose dust away first. Use lens tissue of good quality moistened with a lens cleaner or a small amount of alcohol; then wipe the surface clean with a lens tissue. Solvents such as xylene should NOT be used as cleaner.

3. Do not disassemble objective lenses!

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Bresser GmbH, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

EN 61326
EN 61000-3-2:2005

Product description: Metallurgical reflection-type microscope
Model: BRESSER Science MTL-201

Rhede, 20. July 2007

Bresser GmbH

Helmut Ebbert
Managing director

WARRANTY

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Devices which become defective during the warranty period can be returned to the dealer where the device was bought. The repaired device or a new one will then returned to you. In the case of defects which occur after the end of the warranty period, the devices can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

Important:

Make sure to return the device carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This warranty does not imply any restriction of your statutory rights.

Your dealer: _____

Name: _____

Postcode / City: _____

Street: _____

Telephone: _____

Date of purchase: _____

Signature: _____

Bresser GmbH

Gutenbergstr. 2 · DE-46414 Rhede · Germany

Tel. +49 (0) 28 72 - 80 74-210

Fax +49 (0) 28 72 - 80 74-222

www.bresser.de · service@bresser.de